## LipoClear

REF LC10, LC15, LC40, LCP1

### Cautions

 *Universal Precautions should be followed on all specimens, regardless of whether a specimen is known to contain an infectious agent.*

*If wrong sample volume is used, incorrect results may be obtained.*

### Intended Use

For *in vitro* diagnostic use.

LipoClear is a non-toxic, non-ionic solution used to remove lipoproteins in grossly lipemic specimens, yielding a clear sample suitable for analysis.

### Specimen Requirement

- Fresh serum is the preferred specimen.
- Allow frozen or refrigerated specimens to come to room temperature before processing.
- Analyze cleared specimen within one hour to prevent additional loss of serum proteins.

### Limitations

- LipoClear CANNOT be used to clarify specimens requiring lipid testing.
- LipoClear CANNOT be used to estimate gamma globulin.
- LipoClear CANNOT be used to estimate immunoglobulins.

Specimens treated with LipoClear lose approximately 0.4 - 0.8 g/dl protein, due to the loss of lipoproteins.

Phosphorus was spuriously elevated by about 10% on several analytical systems.

Electrolytes measured on undiluted samples should not be multiplied by appropriate correction factor for each analyte.

### Procedure

- Check product number before proceeding.
- Allow reagent to come to room temperature (let stand for 10 minutes).
- If using LCP1, invert bottle 5 times.
- Proceed with sample addition using the above information for "designated sample volume".

- Shake well or vortex to mix.
- Allow sample/reagent mixture to stand for 5 minutes, then centrifuge.

### Centrifugation

- StatSpin MP, RP, VT: spin on "LipoClear 95seconds" or the "Hard" cycle.
- Express, Express 2 or VET Express: spin on the "Hard" cycle.
- Conventional centrifuges: 2000 x g, for 20 minutes.

Remove cleared sample for analysis.

- If lipids do not form a button at the bottom of the tube, a floating "pedicle" may result.

Pipette past floating pedicle to obtain cleared specimen.

Analyze cleared specimen using laboratory's standard method.

**All results must be multiplied by 1.2 to compensate for the dilution of the original sample.**

- If lipids do not form a button at the bottom of the tube, a floating "pedicle" may result.

Pipette past floating pedicle to obtain cleared specimen.

Analyze cleared specimen using laboratory's standard method.

**All results must be multiplied by 1.2 to compensate for the dilution of the original sample.**

If lipids do not form a button at the bottom of the tube, a floating "pedicle" may result.

Pipette past floating pedicle to obtain cleared specimen.

Analyze cleared specimen using laboratory's standard method.

**All results must be multiplied by 1.2 to compensate for the dilution of the original sample.**

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