


ThermoBrite Exploded View

ThermoBrite



ITEM	PART NUMBER	DESCRIPTION
1	01-004915-001	FUSE 3.0A (120V unit)
	01-004915-002	FUSE 1.6A (240V unit)
2	01-004916-001	FUSE DRAWER (120V unit)
	01-004916-002	FUSE DRAWER (240V unit)
3	01-004917-001	FEET
4	01-004918-001	CABLE 1
5	01-004919-001	CABLE 2
6	01-004920-001	CABLE 3
7	01-004921-001	AC INPUT
8	01-004923-001	ASSY POWER SUPPLY
9	01-004924-001	PCB INTERFACE
10	01-004925-001	ASSY DISPLAY
11	01-004934-001	ASSY HOUSING
12	01-004927-001	ASSY PLENUM & FAN
13	01-004928-001	ASSY SEPARATOR
14	01-004929-001	HINGES
15	01-004935-001	ASSY COVER
16	01-004931-001	COVER GASKET
17	01-006096-001	PLENUM GASKET
Items not shown		
18	01-002532-001	LINE CORD (120)
19	01-004987-001	RS232 MANUAL
20	01-004861-001	RS232 KIT
21	01-004988-001	OPERATORS .MANUAL
22	01-004997-001	BOX WITH FOAM SET
23	01-004983-001	FOAM SET

Part Numbers and availability are subject to change without notice. Contact Customer Service for current information.